



Strengthen Your Culture of Safety with Robust and Actionable Data

PURELL SMARTLINK™ Activity Monitoring System

Perform group-level monitoring with 24/7 data collection, using room sensors that count all activity in and out of patient areas. Data includes healthcare personnel, doctors, patients, and visitors to help form big-picture analysis.

- Automatically track soap and sanitizer events and opportunities at entry and exit points in a monitored area
- Capture robust, accurate, and actionable data for immediate feedback and effective interventions
- Access data in easy-to-use, customizable reports that graphically display performance metrics for staff and leadership
- Combine with PURELLSM Clinician-Based Support to provide customized plans and unbiased assessments of your hand hygiene program
- Use with PURELL SMARTLINK™ Service Alerts to help decrease product waste, improve operational efficiency, and ensure product availability
- Helps you meet Leapfrog Hospital Survey requirements



How It Works

The PURELL SMARTLINK™ Activity Monitoring System uses proprietary technology to communicate with our smart, industry-leading dispensers and patient room activity counters. Battery-powered sensors featuring wireless technology reduce installation complexity.



Activity Counter
Tracks entries and exits in patient rooms and monitored areas




PURELL SMARTLINK Dispensers
Capture soap and sanitizer dispenser activations



Network Gateway
Receives captured events and opportunities and sends data through the Microsoft Azure Cloud

PURELL SMARTLINK Software
Converts near real-time data into actionable reports and visual displays



Free 60-Day Trial

Our unique trial process allows you to validate the AMS technology through interaction with the hardware and software, and experience the expertise of our clinicians through targeted hand hygiene initiatives.

To learn more, call **1-800-321-9647** or email **SMARTLINK@GOJO.com**

PURELL SMARTLINK™ Activity Monitoring System helps you meet **LEAPFROG HOSPITAL SURVEY** requirements with robust 24/7 data collection of hand hygiene opportunities and performance.

©2022 GOJO Industries, Inc. All rights reserved.
33062 (1/2022) | HCE-LIT-PURAMS



Visit **GOJO.com/SMARTLINK** to learn more or email **SMARTLINK@GOJO.com**